





iCHIROPRO – A REVOLUTION FROM BIEN-AIR

This new device, controlled from an iPad and equipped with some of the best performing instruments on the market, paves the way for a new generation of medical appliances. Its innovative solutions make practitioners' lives easier – that's the Bien-Air promise.

YOUR NEW PHYSIODISPENSER

MX-i LED micromotor with maintenance-free ceramic ball bearings which are lubricated for life. Its high torque - even at low speeds - offers unrivalled working comfort.

The **CA 20:1 L Micro-Series internal irrigation contra-angle**, the solution for smooth and efficient operations, is fitted with one of the smallest heads ever designed and a double LED system for uniform lighting regardless of the speed of rotation.

iPad controlling your implantology system offering unprecedented practical functions and guaranteeing unrivalled ease of use.



THE iCHIROPRO iPad APPLICATION

Simplified operations thanks to an intuitive and ergonomically designed application. Available from the Apple App Store, the iChiropro application opens up numerous avenues in terms of customisation options. New innovative functions will regularly be added to further enhance the performance levels of your system.



Multiple users with storage of individual settings



Integrated patient file on the iPad control interface



Complete pre-programmed operating sequences of the largest implant manufacturers



Operation data stored in the form of graphics and tables, attached to the patient file



Save and export operation parameters whilst documenting them with the patient data



iCHIROPRO BASIC SET WITH CA 20:1 L MICRO-SERIES CONTRA-ANGLE	ref. 1700440-001
1 iChiropro basic set	ref. 1700439-001
1 CA 20:1 L Micro-Series contra-angle	ref. 1600692-001

iCHIROPRO BASIC SET	ref. 1700439-001
1 iChiropro control unit	ref. 1600784-001
1 MX-i LED micromotor	ref. 1600755-001
1 MX-i LED micromotor cable	ref. 1600606-001
1 3-button foot pedal	ref. 1600631-001
10 sterile irrigation lines	ref. 1500984-010

www.ichiropro.com